



DarPharma
WE IMPROVE LIVES !

1039

1378

☐ **Instulation Report**
☒ **Service Report**

CSR No.001

Customer name

Alorf Hospital

Date 23/7/2024

Contact person

2210 3545 1310

Date of Purchase:

Warrantee Validity:

Nature of call

Problem reported

Warrantee Visit

Service contract Visit

Device Not Working?

☐ Yes

☒ No

Equipment type:

Sudoscan 2

Serial number:

Call reported by

Date | time

23/7/2024

Service Details

Service rendered

calibration
quality check

100% scan

Engineer's Remarks

the device is working fine

Status after service

☒ Complete

☐ Incomplete

☐ Pending for spares

☐ Under Observation

Customer Feedback

Remarks

Please rate this call by choosing an option

☐ Extremely Satisfied

☐ Satisfied

☐ Dissatisfied

Name

Amr Bawa

Designation

Adel wazeh

Phone | Fax

65991613

Email

Signature

[Signature]

Date

Place



DarPharma
WE IMPROVE LIVES !

ADDRESS: KUWAIT, SALMYIA, BLOCK 3, LAILA COMPLEX, 1ST FLOOR, OFFICE 4

☎ 99004991

☎ 25713133

✉ INFO@DAR-PHARMA.COM

🌐 WWW.DAR-PHARMA.COM



TEST REPORT

23-07-2024

Supervisor_Physician

MODULE ID : 1229 - TYPE : B - FW VER : 519

SOFT VERSION : v3.3.6024 - DB SERVICE DLL VER : v3.3.6024 - KEYBOARD DLL VER : v3.3.6024

EZSCAN DLL VER : 8963

ELECTRODES : HANDS FW VER : v3.0.4532 - FEET FW VER : v3.0.4532

Operator name (if different from physician name)

Electrode LED test



Reference voltage generation & Reference voltage generator loopback linearity

VR1 1	630	503	446	389	332	275	218	98
VR1 2	646	516	458	399	341	282	224	100

Output channel selection & Input channels low noise front and qualification

	0	1	2	3	4	5
	834	834	834	835	835	834
Other channel	Pass	Pass	Pass	Pass	Pass	Pass

Internal resistances calibration

	0	1
	496	516

Channel sensitivity control

RE1	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
	58	82	107	136	166	194	229	270	316	357	407	458	517	576	641	700
	-120								22							

Input channel linearity control (%) & Input channel absolute linearity control

%	0	10	20	30	40	50	60	70	80	90	100
	0	105	209	315	400	516	616	709	814	912	1023

Electrodes shortcircuit test

%	0	0	98	98
	Pass	Pass	Pass	Pass

Run test

Report

Manual test

More details

?