

ATCO Technology Co.

Tel.: (+965) 24743017 / 24567725
Fax: (+965) 24743019 / 24567792
P.O. Box: 2355 Jahra, 01025 - Kuwait
Email: sales@atco-kuwait.com
www.atco-kuwait.com
CR NO . 336218



شركة أتكو التكنولوجية

تلفون : ٢٤٧٤٣٠١٧ / ٢٤٥٦٧٧٢٥ (+٩٦٥)
فاكس : ٢٤٧٤٣٠١٩ / ٢٤٥٦٧٧٩٢ (+٩٦٥)
ص.ب : ٢٣٥٥ الجهراء, ٠١٠٢٥ - الكويت
Email: sales@atco-kuwait.com
www.atco-kuwait.com
س.ت. ٣٣٦٢١٨

SERVICE / INSTALLATION REPORT

No.: 55711

Customer	Aloof Hospital	Equipment	Endoscopy washer
		Model No.	SERIE 3
		Serial No.	25186
Authorized Person		Manufacturer	Soluscope
Customer Complaint		Location	Endoscopy floor 2
Error - 111 - 31		Contract / Order No.	
Services / Check Carried out		☐ Installation	☐ Warranty
☐ Inspection / Cleaning	☐ Optics Checks	☒ Repair	☐ On Charge
☐ Calibration / Performance	☐ Pneumatic Checks	☐ Maintenance	☐ Contract
☐ Interface / Communication	☐ Hydraulic Checks	☐ Upgrade	☐ _____
☐ Safety Checks	☐ Operational Error		
☐ Electrical / Electronics	☐ Mechanical Checks		
☐ Accessories Check	☐ _____		
		Time Spent	

Remarks Service call attendance.

We inspect the unit and found that while performing the cycle error 111 is occur. It is due to the air flow in dosing chamber. We remove the dosing chamber connection and release the air and reconnect. We tried one cycle and it successfully done. Checked and confirmed that the unit is working properly.

No.	Spare Parts Used	Part No.	Qty.
	While changing the chemical bottle		Kindly
	close the bottle lid properly.		

Service Engineer

ATCO Technology Co.

Signature

Date 29/05/2024

Authorized Person

Signature

Date

Customer / User

Signature

Date