

[illegible]

- * Preventive Maintenance (PM)
- * Checked the Instrument
- * Checked and cleaned the Internal parts
- * Cleaned Transducer by brush and RBC clogs Removed
- * Flush Reaction Chambers with Cellclean Manually.
- * Drained Reaction Chambers and Drained Waste Fluid Chamber.
- * Perform Rinse Flowcell and Air Bubbles Removed.
- * Checked FCM Syringe, Sensors and Motors.

REMARKS	
* PM Done.	
JOB COMPLETED ☒ Yes ☐ No	AWAITING PARTS ☐ Yes ☒ No

ATC ENGINEER SIGNATURE <i>Euseb</i> Fasa Diab.	USER SIGNATURE Dr. Arafat Samy	AUTHORIZED SIGNATURE & STAMP دكتور ايهاب محمد عبد الحليم 22020120560
---	--	---