

Atco



MEDICAL EQUIPMENT REPAIR REQUEST

No. _____

Date: 10/8/24

Requested by Rage

Department GA Location Ward

SR. No.	Machine	Model	Complaint
<u>0770320120028</u>	<u>DVT machine</u>	<u>DVT 700</u>	<u>Not working while touching on display.</u>

Check within

 Week

ASAP

Engineer Comments

End Comments

Requested by

Checked By

Approved by

ATCO Technology Co.

Tel.: (+965) 24743917 / 24567725

Fax: (+965) 9 / 24567792

P.O. Box: 201025 - Kuwait

Email: sales@atco-kuwait.com

www.atco-kuwait.com

CR NO . 336218



شركة أتكو التكنولوجية

تلفون : ٢٤٧٤٣٠١٧ / ٢٤٥٦٧٧٢٥ (+٩٦٥)

فاكس : ٢٤٧٤٣٠١٩ / ٢٤٥٦٧٧٩٢ (+٩٦٥)

ص.ب : ٢٣٥٥ الجبراء, ٠١٠٢٥ - الكويت

Email: sales@atco-kuwait.com

www.atco-kuwait.com

س.ت. ٢٣٦٢٨

SERVICE / INSTALLATION REPORT

No.: 61994

Customer	Al orf Hospital	Equipment	DVT Machine
		Model No.	7700
		Serial No.	D770B20120024
Authorized Person		Manufacturer	Wojin
Customer Complaint		Location	CEA WARD
Touch panel not working		Contract / Order No.	
Services / Check Carried out ☐ Inspection / Cleaning ☐ Calibration / Performance ☐ Interface / Communication ☐ Safety Checks ☐ Electrical / Electronics ☐ Accessories Check ☐ Optics Checks ☐ Pneumatic Checks ☐ Hydraulic Checks ☐ Operational Error ☐ Mechanical Checks ☐ _____		☐ Installation ☐ Warranty ☒ Repair ☐ On Charge ☐ Maintenance ☐ Contract ☐ Upgrade ☐ _____	
		Time Spent	

Remarks

Service call attended.
 we inspect the unit, we found that the touch panel display not working
 we replaced the touch panel with from the spare unit in stock.
 we tested the connection and keys and we checked the unit is working properly.

No.	Spare Parts Used	Part No.	Qty.

Service Engineer Nixon / Mousa ATCO Technology Co. Signature Date 28-08-2024	Authorized Person Signature Date	Customer / User Signature Date
--	---	---