

ATCO Technology Co.

Tel.: (+965) 24743017 / 24567725

Fax: (+965) 24743019 / 24567792

P.O. Box: 2355 Jahra, 01025 - Kuwait

Email: sales@atco-kuwait.com

www.atco-kuwait.com

CR NO . 336218



log no 2591

شركة اتكو التكنولوجية

تلفون : ٢٤٧٤٣٠١٧ / ٢٤٥٦٧٧٢٥ (+٩٦٥)

فاكس : ٢٤٧٤٣٠١٩ / ٢٤٥٦٧٧٩٢ (+٩٦٥)

ص.ب : ٢٣٥٥ الجهراء ٠١٠٢٥ - الكويت

Email: sales@atco-kuwait.com

www.atco-kuwait.com

س.ت. ٣٣٦٢١٨

SERVICE / INSTALLATION REPORT

No.: 60984

Customer	Alor Hospital	Equipment	Autoclave
		Model No.	Selectomat PL
		Serial No.	220823
Authorized Person		Manufacturer	Mmm
Customer Complaint		Location	CSSD
Error	cooling water temp high	Contract / Order No.	
Services / Check Carried out		☐ Installation	☐ Warranty
☐ Inspection / Cleaning	☐ Optics Checks	☒ Repair	☐ On Charge
☐ Calibration / Performance	☐ Pneumatic Checks	☐ Maintenance	☐ Contract
☐ Interface / Communication	☐ Hydraulic Checks	☐ Upgrade	☐ _____
☐ Safety Checks	☐ Operational Error	Time Spent	
☐ Electrical / Electronics	☐ Mechanical Checks		
☐ Accessories Check	☐ _____		

Remarks Service call attended. we inspect the unit and found that there is no cool water supply to the machine. we checked the chiller side tank and found that one motor is tripped off. we restart the motor and checked the water flow, now the cool water supply is ok. we performed one cycle and found that the unit is working properly.

No.	Spare Parts Used	Part No.	Qty.
	NB:- Informed maintenance team to check the motor and chiller regularly for the better performance of machine.		

Service Engineer	Authorized Person	Customer / User
Nidion ATCO Technology Co.		
Signature	Signature	Signature
Date 15/07/2024	Date	Date