

ATCO Technology Co.

Tel.: (+965) 24743017 / 24567725

Fax: (+965) 24743019 / 24567792

P.O. Box: 2355 Jahra, 01025 - Kuwait

Email: sales@atco-kuwait.com

www.atco-kuwait.com

CR NO . 336218



شركة أتكوتكنولوجيه

تلفون : ٢٤٧٧٧٢٥ / ٢٤٧٤٣٠١٧ (+٩٦٥)

فاكس : ٢٤٥٦٧٧٩٢ / ٢٤٧٤٣٠١٩ (+٩٦٥)

ص.ب : ٢٣٥٥ الجبراء , ٠١٠٢٥ - الكويت

Email: sales@atco-kuwait.com

www.atco-kuwait.com

س.ت. ٢٣٦٢١٨

SERVICE / INSTALLATION REPORT

No.: 63786

Customer	Aloof Hospital	Equipment	patient monitor-
		Model No.	Elite V6, XMI module
		Serial No.	260801-M1971/220001
Authorized Person		Manufacturer	Edan
Customer Complaint		Location	OR
	MBP not working	Contract / Order No.	
Services / Check Carried out		☐ Installation	☐ Warranty
☐ Inspection / Cleaning	☐ Optics Checks	☒ Repair	☐ On Charge
☐ Calibration / Performance	☐ Pneumatic Checks	☐ Maintenance	☐ Contract
☐ Interface / Communication	☐ Hydraulic Checks	☐ Upgrade	☐ _____
☐ Safety Checks	☐ Operational Error	Time Spent	
☐ Electrical / Electronics	☐ Mechanical Checks		
☐ Accessories Check	☐ _____		

Remarks Service call attended
we inspect the unit and checked the MBP in
xmi module. MBP is not going up, it blocking in middle
we swap the other module and its working properly.
So the xmi module is faulty. we will check with
manufacture and update you as soon as possible

No.	Spare Parts Used	Part No.	Qty.
	Need repair		
	checking done		

Service Engineer Mixon ATCO Technology Co. Signature Date 22/12/2024	Authorized Person Signature Date 21/12	Customer / User Signature Date
--	--	--------------------------------------