

ATCO Technology Co.

Tel: (+965) 24743017 / 24567725

Fax: (+965) 24743019 / 24567792

P.O. Box: 2355 Jahra, 01025 - Kuwait

Email: sales@atco-kuwait.com

www.atco-kuwait.com

CR NO. 336218



logno-2592, 2591

شركة اتكو التكنولوجية

تلفون: ٢٤٧٤٣٠١٧ / ٢٤٥٦٧٧٢٥ (+٩٦٥)

فاكس: ٢٤٧٤٣٠١٩ / ٢٤٥٦٧٧٩٢ (+٩٦٥)

ص.ب: ٢٣٥٥ الجهرام, ٠١٠٢٥ - الكويت

Email: sales@atco-kuwait.com

www.atco-kuwait.com

س.ت. ٣٣٦٢١٨

SERVICE / INSTALLATION REPORT

No.: 58334

Customer	Alorf Hospital		
Equipment	CSSD Autoclave		
Model No.	Seltomat PL		
Serial No.	220824 & 220823		
Authorized Person	Manufacturer M M M		
Customer Complaint	Location CSSD		
Error:- coding water temp high.	Contract / Order No.		
Services / Check Carried out	☐ Installation ☐ Warranty		
☐ Inspection / Cleaning	☒ Repair ☐ On Charge		
☐ Calibration / Performance	☐ Maintenance ☐ Contract		
☐ Interface / Communication	☐ Upgrade ☐		
☐ Safety Checks	Time Spent		
☐ Electrical / Electronics			
☐ Accessories Check			
Remarks	Service call attended.		
when we inspect the unit the cycle before run having the error cooling water temp high.			
we checked the input water supply and its temp its all are in the limit - It seems that at the time of cycle some one the water supply is interrupted.			
we tried new cycle on both unit.			
No.	Spare Parts Used	Part No.	Qty.
	we found that the unit is working properly		
under Warranty			
Service Engineer		Authorized Person	
Nixon		22/2/24	
ATCO Technology Co.			
Signature		Signature	
Date 21/02/2024		Date	
		Customer / User	
		Signature	
		Date	