

Customer Report Form



Case number CAS-0021210349

Order number WO-29829014 | Visit Date 09/05/2024

Instrument Details

Instrument/Module cobas b 221<6>=Roche OMNI
S6 system

Serial Number 19239

Additional Details

PO

Customer Details

Customer Number 0050250925

Customer Name Dr. Mansour Alorf Hospital

Street Address Mansour Saad Al Bannaq
Street

Zip 01025 - Jahra

Contact Name Kristina Marcelino

Contact Phone +96565991613

Service Activity Code Service Visit

Purpose of Visit

CUSTOMER ALLEGATION

=====

Created by: Gina Bisenio, biomed@alorf.com (+96565991613)

Preferred callback time:

Instrument: b 221 6 Roche OMNI S6 system (Blood Gas/Electrolytes) (Serial: 19239)

Alarm code: hw test error

Issue type: OTHER

Lot no/Vendor: not applicable

Initial description: hw test error

Issue timeline:

1: Issue created (05.09.2024 15:55 GMT+00:00)

ACTIONS

=====

Type your descriptive analysis here

CONCLUSION

=====

Log a call or note

Performed Activities

The customer not using ISE and it was de activated.

Customer re activated it causing temp alarm. so the ISE parameters were deactivated again and machine rebooted and working fine.

Time Report

Category	Start	End	Hours
05913616001-Service labour time	09/05/2024 11:30 AM	09/05/2024 2:00 PM	2.50

Signature or Stamp

Customer
Gina Bisenio

Roche Support
Sunash Shukoor

Sign Here

WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO
WO-20820014 12/10/2024, 3:29 PMWO

By signing this document

- Roche Support person confirms the above mentioned work and checks have been completed..
- Customer representative confirms they have reviewed the work carried out as described and accepts the system for use.
- Thank you for your trust.

Address: Al Mirqab, Al Soor Street, Block 1, Burj Al Shaya, Plot 6, Floor 5, Office 11.

P.O. Box: 467, Safat - 13005 Kuwait

Hotline: 25770360

Fax: 22059262

Customer Assistance: