

Customer Alorif Hospital		Equipment CSSD Autoclave	
		Model No. Selectomat PL	
		Serial No. 220823	
Authorized Person		Manufacturer MMM	
Customer Complaint		Location CSSD	
Error:- cooling water temperature High		Contract / Order No.	
		☐ Installation ☒ Repair ☐ Maintenance ☐ Upgrade ☐ Warranty ☐ On Charge ☐ Contract ☐ _____	
Services / Check Carried out ☐ Inspection / Cleaning ☐ Calibration / Performance ☐ Interface / Communication ☐ Safety Checks ☐ Electrical / Electronics ☐ Accessories Check ☐ Optics Checks ☐ Pneumatic Checks ☐ Hydraulic Checks ☐ Operational Error ☐ Mechanical Checks ☐ _____		Time Spent	
Remarks Service call attended. we inspect the unit and found that the cooling water temperature is too high. we checked the Chiller unit and found that the unit is not working. Without chiller the two autoclave not working properly. we informed maintenance team to fix the issue as soon as possible.			
No.	Spare Parts Used	Part No.	Qty.
Service Engineer Midcom ATCO Technology Co. Signature Date 08/07/2024		Authorized Person Signature Date 16/9/24	
		Customer / User Signature Date	